



RANDOLPH AUSTIN COMPANY

Acceptable Sterilization Methods – Cilran™ Tubing

Steam Sterilization

- The following times can be suggested for the maximum autoclave cycles for Cilran™ Tubing:
 - 250° F (121° C) for 17 minutes
 - 240° F (116° C) for 27 minutes.
- There may be some dimensional change due to the elevated temperatures of autoclaving. This is because of the stress relieving that may occur. It is strongly recommended that there be no outside mechanical stress on the tubing while it is being autoclaved.

Gamma Sterilization

- Cilran™ can be sterilized using gamma radiation without suffering a large loss in the physical properties, even after extended exposure.

Product Disclaimer

- In regards to Cilran™ tubing, Randolph Austin Company makes no warranties of any kind, either expressed or implied, including but not limited to warranties of merchantability, fitness for a particular purpose, of title, or of non-infringement of third party rights. Use of the product by a user as is at user's risk.